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# **Standard Operational Procedure for External Quality Assessment in TB microscopy**

## Contents

	Page
<b>Foreword .....</b>	<b>3</b>
<b>Acknowledgement .....</b>	<b>4</b>
<b>1. Introduction .....</b>	<b>5</b>
<b>2. Glossary of Terms:.....</b>	<b>6</b>
<b>3. EQA Method applied in Camodial.....</b>	<b>9</b>
<b>3.1 TB Laboratory network.....</b>	<b>9</b>
<b>3.2 EQA network .....</b>	<b>10</b>
<b>3.3 Flowchart of EQA activities.....</b>	<b>11</b>
<b>3.4 Blinded rechecking .....</b>	<b>12</b>
<b>3.5 On site evaluation .....</b>	<b>13</b>
<b>3.6 EQA quarterly work shop.....</b>	<b>15</b>
<b>3.7 Reporting to CENAT.....</b>	<b>16</b>
<b>4. Reference.....</b>	<b>17</b>
<b>5. Appendix.....</b>	<b>18</b>

## Foreword

The Cambodia National Tuberculosis Program (NTP) has been established to particularly detect infectious cases of pulmonary tuberculosis based on the bacteriological examination, to monitor treatment progress and document cure at the end of treatment by means of sputum smear microscopy. Bacteriological services also contribute to the diagnosis of pulmonary and extra-pulmonary tuberculosis by appropriate methods. Therefore, standardization of the basic tuberculosis bacteriological techniques becomes a mandatory requirement.

A reliable laboratory service is one that is cost-effective and provides results that are consistently accurate. These demands can be met only through commitment to quality assurance. A key component of quality assurance for sputum smear microscopy services is External Quality Assessment (EQA), which is a process to assess and monitor laboratory performance.

In Cambodia, for EQA the NTP focuses on two main methods namely blinded rechecking and on site evaluation by supervisory visit.

I hoped that this EQA SOP will be implemented in Cambodia as a means of ensuring the high quality of the NTP. In addition, I believe that this SOP is also of great importance for National Reference Tuberculosis Laboratory including Provincial TB Laboratory Supervisors, OD supervisors and microscopists to follow.

Phnom Penh, 29 April, 2010

**National Tuberculosis and Leprosy Control  
Director**

Dr. Mao Tan Eang

## **Acknowledgement**

On behalf of the working group for the formulation of the EQA SOP (External Quality Assessment Standard Operational Procedure) in Cambodia, I would like to express our sincere gratitude to all members of the group who have participated and made the EQA SOP successful.

We wish to particularly thank Ms. Hiroko Matsumoto of JICA-National TB Control Project for her great effort in leading, consulting with all NTP laboratory staffs (especially, EQA staffs and provincial TB laboratory supervisors), and creating this fundamental tool.

We would also like to emphasize that without strong support, full and active participation of JICA TBV Control Project, NTP, and all partners' staffs involved, the finalization of this EQA SOP document could not successfully take place.

## 1. Introduction

Since the diagnosis of tuberculosis relies on sputum smear microscopy, maintenance of quality of microscopy is vital to ensure reliable and accurate microscopy service of NTP. EQA was started to continuously improve the reliable, efficient and sustainable of TB laboratory services. EQA includes (1) on-site evaluation of the laboratory to review QC and should include on-site rereading of smears; (2) blinded rechecking ; and (3) panel testing . EQA also allows participant laboratories to assess their capabilities by comparing their results with those obtained in other laboratories in the network (intermediate and central laboratory) through panel testing and rechecking. In Cambodia, blinded rechecking and on site evaluation were chosen for the main methods for EQA implementation.

The purposes of EQA are the followings;

- To improve quality of services
- To identify cause of errors
- To maintain networking between the laboratories
- To compare laboratory results to intermediate and central level
- To Share experience including good and bad

This SOP is prepared for those who involved in EQA activities. i.e. Assessor, Controller, Provincial TB laboratory supervisor (PLS), OD supervisor and Microscopist.

This EQA SOP is prepared for blind rechecking, on site evaluation, EQA work shop in Cambodia.

## 2. Glossary of Terms

**Quality Assurance (QA):** System designed to continuously improve the reliability and efficiency of laboratory services. It includes internal quality control, external quality assessment, and quality improvement.<sup>1)</sup>

**Quality Control (QC):** Also called Internal Quality Assurance, includes all means by which the TB smear microscopy laboratory controls operation, including instrument checks and checking new lots of staining solutions.<sup>1)</sup>

**External Quality Assessment (EQA):** A process which allows participant laboratories to assess their capabilities by comparing their results with those in other laboratories in the network (intermediate and central laboratory) through panel testing and blinded rechecking. EQA also includes on-site evaluation of the laboratory to review quality of performance and should include on-site rereading of smears. EQA is an expansion of the proficiency testing as described by IUATLD.<sup>1)</sup>

**Quality Improvement (QI):** A process by which the components of smear microscopy diagnostic services are analyzed with the aim of looking for ways to permanently remove obstacles to success. Data collection, data analysis, and creative problem solving are the key components of this process. It involves continued monitoring, identifying defects, followed by remedial action including retraining when needed, to prevent recurrence of problems. QI often relies on effective on-site evaluation visits.<sup>1)</sup>

**Panel Testing:** Sending stained and/or unstained smears from the reference laboratory to the peripheral or intermediate laboratory to check proficiency in reading and reporting.

Panel testing is equivalent to the WHO definition of proficiency testing. The term panel testing is used in these guidelines in order to eliminate the confusion over the different definitions of proficiency testing.<sup>1)</sup>

**Rechecking:** Sending smears from the peripheral laboratory to a reference laboratory (intermediate or central laboratory) for rereading. These guidelines recommend that rechecking is always blinded, ensuring that the controller does not know the results from the peripheral laboratory. In other documents, this may also be referred to as rereading.<sup>1)</sup>

**Blinded Rechecking:** Sending smears from the peripheral laboratory to a reference laboratory (intermediate or central laboratory) for rereading by blinded, ensuring that the controller does not know the results from the peripheral laboratory.

**Assessor:** Term used to describe the supervisor of laboratories. Assessor is responsible for slide collection and selection, umpire reading for the discordant slide between tested laboratory and controllers, on site evaluation, organizing EQA workshop and report the EQA results to CENAT (national TB reference laboratory)

**Controller:** Term used to describe the technician responsible for rechecking slides.<sup>1)</sup>

**Major error:** This type of error is considered the most critical since it has the highest potential impact on patient management, and can result in an incorrect diagnosis or improper management of a patient. Major errors may indicate gross technical deficiencies, and include both High False Positive and High False Negative errors.<sup>1)</sup>

**High False Positive (HFP):** A negative smear that is misread as 1+ to 3+ positive. This is a major error.<sup>1)</sup>

**High False Negative (HFN):** A 1+ to 3+ positive smear that is misread as negative. This is a major error.<sup>1)</sup>

**Minor error:** In clinical practice, these errors may have some impact on patient management. However, for the purpose of evaluating laboratory performance, this type of error is considered less serious, because of inherent limitations in consistently detecting a few AFB that may be unequally distributed within a smear. The frequency of minor errors may indicate technical deficiencies.<sup>1)</sup> *In Cambodia LFP and LFN are minor error and QE is not include.*

**Low False Positive (LFP):** Previously called a scanty false positive.<sup>1)</sup> A negative smear is misread as a low (1-29 AFB/300 fields) positive. This type of minor error occurs occasionally even in laboratories that are performing well.

**Low False Negative (LFN):** Previously called a scanty false negative.<sup>1)</sup> A low (1-29 AFB/100 fields) positive smear that is misread as negative. This type of minor error occurs occasionally even in laboratories that are performing well.

**Low Positive:** Term used in this document to describe 1-9 acid-fast bacilli per 100 fields, which is the WHO/IUATLD standard for quantitation. (1-29 AFB/300 VF in Cambodia) These results are reported to the physician as exact number of AFB seen. It is up to the physician and the NTP to decide if this represents a case or not. It previously referred to as a scanty positive.<sup>1)</sup>

**Feedback:** Process of communicating results of EQA to the original laboratory, including suggestions for possible causes of errors and remedies.<sup>1)</sup>

**Unaccepted laboratory:** Laboratory in which major error or more than one minor error is found during a period of the time (ordinary one quarter)

### 3. EQA Methods Applied in Cambodia

#### 3.1 TB Laboratory network

It is important to provide TB smear microscopy services that are accessible to the entire population, yet to maintain an acceptable level of technical proficiency. To accomplish this objective, a network of laboratories with competency in acid-fast sputum smear microscopy, supported by provincial TB laboratory supervisor, and overseen by a National TB reference laboratory, is required.


Figure 1. TB laboratory network

### 3.2 EQA network



Figure 2. EQA network

### 3.3 Flowchart of EQA activities

### 3.4 Blinded Rechecking

Blinded rechecking of a sample of routine smears from the peripheral sites and intermediate labs by controllers at a higher level laboratory is considered the best method for evaluation performance and providing motivation to staff for improvement.

#### 1. Slide storage

All examined slides should be store serially in the slide box until slide selection will be done

#### 2. Slide selection


- (1) Twenty slides per laboratory are selected for blind-rechecking in each quarter
- (2) Count the total slide number in the test quarter for blind rechecking and quarterly activity report. If total slide number is less than 20, collect all slides.
- (3) Calculate the sequence number (total slide number divide by 20) and fill it in the slide selection sheet (Appendix 1)
- (4) Chose the starting number randomly and it should be less than sequence number.
- (5) Select the slides from the register book but not from slide box.
- (6) Fill the slide number and results to the slide selection sheet (Appendix 1)
- (7) Chose the slides according to the slide selection sheet. In the case of needed slide is missing or broken, chose the next slide and erase, using pen, the original selection and put the replace slide number and result. Note the reason in the remark column
- (8) Verify the slide number and its results in the slide selection sheet. After confirmation with Microscopist, obtain signature from him/her ]

3. Cross checking
  - (1) Assessor transfer only slide number to feedback sheet (Appendix 2)
  - (2) Assessor give the feedback sheet and selected slide to controller without letting him/her know the results (blinded)
  - (3) Controller re-stains all slides by Z-N staining method (To avoid fading).
  - (4) Controller fills the feed-back sheet after reading the slides: date of examination, results, and smear assessment such as specimen quality, thickness, size and evenness, by ticking in each category.
4. Analysis of results (Assessor)
  - (1) Check the feed-back sheet which is filled by controller. If it is not complete, ask controller to complete it.
  - (2) Fill the microscopist's results to the feed-back sheet
  - (3) Compare controller's results and microscopist's (original results)
  - (4) Re-check the slides if discordant slides are found.
  - (5) Calculate the percentage of each smear assessment points
  - (6) Fill correlation table, according to the results of microscopist and controller/assessor (Appendix 5)
  - (7) Fill error classification table according to correlation table
  - (8) Fill comments as finding of the results
  - (9) Date the second cross-check and date the last analysis with signature.

### 3.5 On Site Evaluation

Visits to the peripheral laboratories by trained laboratory personnel from the reference or intermediate laboratory are essential if performance is to be improved or maintained at a high standard. These visits allow for the observation of worker performance under actual conditions, including condition of equipment, laboratory safety, adequacy of supplies, and the process for smearing, staining, reading, recording and reporting. Stained smears can be reviewed during the visit. When problems are detected, solutions can be suggested and potentially implemented immediately.

1. Make a plan to visit laboratories which have experienced major error or more than one minor errors
2. According to “Criteria of on site evaluation” (Figure 4), select the laboratory where Assessor has to perform on site evaluation
3. Required forms and equipments are listed below,
  - slide selection sheet (Appendix 1)
  - selected slides (at least discordant slides)
  - feed-back sheet (Appendix 2)
  - blank on site evaluation sheet (Appendix3)
  - the latest visit of on site evaluation sheet (Appendix 3)

**Figure 4.** Criteria of on site evaluation

\* In the case of one minor error, its next quarter results should be followed up.

4. Confirm the discordant results in register book
5. Ask microscopist to read discordant slides together with assessor.
6. Examine at least 5 positive slides at the place.
7. Observe condition and practice
8. Find the problem, cause of problems and also good points
9. Make recommendation to laboratory, obtain signature of microscopist and assessor on the on site evaluation sheet, and leave the carbon copy of on site evaluation sheet (Appendix 3)
10. Summary Report to head of institution
11. Follow up should be done in necessary timing

### 3.6 EQA Quarterly Workshop

EQA quarterly work shop is formulated to share information between stake-holders. There is some possibility to have errors even the laboratory which has no error in EQA samples. This is the opportunity to discuss cause or risk of error between all laboratory workers in the province.

1. Summarize all laboratory results in a province. Fill quarterly laboratory activity and results of blinded check in EQA report sheet (Appendix 4)
  - **Agreement rate:** “correctly reading slides” divide by “total slide number” and multiply 100. In case of Cambodia, “correctly reading slide” means correctly read as positive or negative by qualitative and it is not for quantitative.
  - **Unacceptable laboratory:** The laboratory which found 1 major error or more than 1 minor error called “unacceptable performance laboratory” or “unacceptable laboratory”

Followings are also able to analyze;

- **False negative rate:** In Cambodia, False negative divide by tested laboratory negative multiply 100.
- **False positive rate:** In Cambodia, False positive divide by tested laboratory positive multiply 100
- **Percentage of unacceptable laboratory:** Number of unacceptable laboratory divide by total laboratory number multiply 100

2. Hold EQA quarterly meeting by Provincial Health Department /Or CENAT depend on available budget and facilities.
3. Attendance listed below:
  - Director of Provincial Health Department (PHD)
  - Chief Technical Bureau in PHD
  - CENAT laboratory supervisor
  - Provincial TB medical supervisor
  - Provincial TB laboratory supervisor (reporter)
  - OD supervisors
  - Chief of General laboratory in Provincial hospital
  - Controller
  - Microscopists
  - some other HC staff who involve in TB activities
4. Report EQA results of the last quarter
5. Report the results of on site evaluation and share experiences
6. Discuss about critical problems or findings  
Assessor submit the copy of EQA quarterly workshop to CENAT, national TB reference laboratory

### 3.7 Reporting to CENAT

Submit the copy of reports to CENAT, National TB Reference Laboratory

- Slide selection sheet (Appendix 1)
- Feed back sheets of all laboratories in province (Appendix 2)
- On site evaluation sheets (Appendix 3)
- EQA report sheet (Appendix 4)
- Quarterly laboratory activity report

#### 4. Reference

- 1) IUATLD, WHO, JATA, KNCV, CDC, APHL, External Quality Assessment for AFB Smear Microscopy, APHL, Sept. 2002),
- 2) WHO, Quality Assurance of Sputum Microscopy in DOTS programme, Regional Guideline for countries in the Western Pacific , WHO 2003
- 3) Cambodia EQA Guideline (Draft/ July 2006)

## **5. Appendix**

**Appendix 1 Slide selection sheet**

**Appendix 2 EQA Feedback sheet**

**Appendix 3 On site evaluation sheet**

**Appendix 4 EQA report sheet**

**Appendix 5 Correlation table and classification of error**

**Appendix 6 Definition and Criteria for acceptable smear slide preparation**

**Appendix 7 Member of EQA SOP working grou**

### Slide Selection Sheet

Date of selection	Quarter	Year
Province	OD	
Name of Lab		
Name of Lab tech		
Name of supervisor		

Diagnose				Follow-up			
	Positive	Negative	Total		Positive	Negative	Total
D1				C2,C3,C4			
D2				C5			
D3				C6,C7,C8			
Total			<sup>a</sup>	Total			<sup>b</sup>

Total # of slide in the quarter :	<sup>c</sup>	= (	<sup>a</sup>	+	<sup>b</sup>	)
Sequence number:		=	<sup>c</sup>	/	20	
Starting number:		randomly and less than sequence #				

**Selected Slides**

#	Slide Number	Result	Remarks
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
11			
12			
13			
14			
15			
16			
17			
18			
19			
20			

Signature of Microscopist

Signature of Supervisor

### EQA Feed-Back Sheet

<b>Date of cross check:</b>					<b>Lab code:</b>											
<b>Province:</b>		<b>OD:</b>			<b>Laboratory:</b>				<b>QA Center:</b>				<b>Q:</b>	<b>Year:</b>		
Slide No	Original Slide No.	Results			Specimen Quality			Thickness			Size			Evenness		
		Lab	C.C	Ass	good	saliv	noth	good	thick	thin	good	big	small	good	bad	
1																
2																
3																
4																
5																
6																
7																
8																
9																
10																
11																
12																
13																
14																
15																
<b>total</b>																
<b>%</b>																

Results of crosschecker/ assessor	Results of tested laboratory						
		0	1-29	1+	2+	3+	Total
	0						
	1-29						
	1+						
	2+						
	3+						
	Total						

Summary of errors	
High False negative	
High False positive	
Minor error	
Total	

**Comments**

**Signature of microscopist**

**Date**

**Signature of assessor**

**On site evaluation sheet**

Name of laboratory:

OD:

Province:

Date of visit :

Date of last visit:

Name of microscopist :

Finding

Recommendation

Follow-up :

Date : .....

Name and signature of microscopist

Name and signature of supervisor

**EQA report** (quarter:                      year:                      )

No	OD	Name of laboratory	True Positive	No of False positive				TRUE Negative	No of False negative			Error classification		Agreement rate	Total Slide check	Total slide examined
				BAR	1+	2+	3+		BAR	1+	2+	3+	minor			
1																
2																
3																
4																
5																
6																
7																
8																
9																
10																
11																
12																
13																
14																
15																
<b>Total</b>																

## Correlation table and Classification of Errors

		Microscopist					
		0	1-29	1+	2+	3+	Total
Controller/ Assessor	0	A	W	X	Y	Z	
	1-29	R	B				
	1+	S		C			
	2+	T			D		
	3+	U				E	
	Total						

Letter located at Correlation Table	Results		Error Classification	
	Lab	CC Ass		
A	0	0	No Errors	Correct
R	0	1-29	LFN (Low False Negative)	Minor Error
S	0	1+	HFN (High False Negative)	Major Error
T	0	2+		
U	0	3+		
W	1-29	0	LFP (Low False Positive)	Minor Error
X	1+	0	HFP (High False Positive)	Major Error
Y	2+	0		
Z	3+	0		

## Definition of Good / Acceptable Quality of Smear Slides Preparation

Stained smear slides can be evaluated whether if they are good or poor in terms of the dominance of the following checkpoints in the smear area macroscopically and microscopically.

- 1) **Specimen Quality:** The presence of dust cell (macrophage) or presence of more than 25 leucocytes per field at total magnification of x100 (objective x 10) are observed. A least 10 fields should be observed.
- 2) **Smear Size:** Approximately 2 x 3 cm in size (1.8-2.0 x 2.8-3.0 cm in size)
- 3) **Smear Thickness:** The whole depth of the smear layer can be focused sharply in each field
- 4) **Evenness:** Surface of smear have to spread evenly by small coiling shape and not to slough off

Note: Percentage of good smear should be equal or more than 80 % for all assessment points.

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