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### ឧទាហរណ៍១: Line-list

No.	Hospital name	Specimen types	Patient ID	Patient information				Results				
				In vs out patients	Age	Gender	Admission date	Date sample collected	Cultures	Pathogens	AST	Resistance mechanism if known
1		Blood	1	Inpatient	30	M	4/12/2017	4/15/2017	Positive	E. coli		
2		Urine	1	Inpatient	30	M	4/12/2017	4/15/2017	Negative	Nothing		
3		Stools	1	Inpatient	30	M	4/12/2017	4/16/2017	Positive	E. coli		
4		Blood	2	Inpatient	62	F	4/7/2017	4/13/2017	Positive	A. baumannii		
5		Urine	2	Inpatient	62	F	4/7/2017	4/13/2017	Negative	Nothing		
6		Blood	3	Inpatient	55	F	4/2/2017	4/8/2017	Positive	E. coli		
7		Blood	4	Inpatient	60	M	4/1/2017	4/5/2017	Negative	Nothing		
8		Vaginal	5	Inpatient	35	F	4/10/2017	4/15/2017	Negative	Nothing		
9		Urine	6	Inpatient	70	F	4/3/2017	4/5/2017	Positive	E. coli		
10		Stools	6	Inpatient	70	F	4/3/2017	4/6/2017	Negative	Nothing		
11		Blood	7	Inpatient	20	M	4/11/2017	4/15/2017	Positive	A. baumannii		
12		Vaginal	8	Outpatient	25	F	NA	4/6/2017	Positive	N. gonorrhoeae		
13		Urine	9	Inpatient	80	F	4/8/2017	4/12/2017	Negative	Nothing		
14		Stools	10	Inpatient	2	M	4/19/2017	4/26/2017	Positive	Shigella sp.		

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**Epidemiology Form (ANNEX 1)**

AMR Sentinel Sites Tier 2:  AHC  BTB  TKO  SCH

Monthly: \_\_\_\_\_ (From: \_\_\_\_\_ / \_\_\_\_\_ / 20 to \_\_\_\_\_ / \_\_\_\_\_ / 20 )

		Number of admitted inpatient to Hospital																		
		All ages		0-28 days		29 days-11 months		1-4 years		5-14 years		15-24 years		25-49 years		50-64 years		≥65 years		
Gender		M	F	M	F	M	F	M	F	M	F	M	F	M	F	M	F	M	F	
<b>IPD</b>																				
<b>Blood</b>																				
<b>Blood positive</b>																				
<b>CSF</b>																				
<b>CSF positive</b>																				

**Hospital/Laboratory Request Form**

SRH  KCH  NPH  CAL

**Hospital Service/ ward:**.....

**Hospital/Laboratory contact:** .....

1. Age:..... BoD: .....Days .....Months .....Years
2. Sex:  Male  Female
3. Patient's hospital ID#:.....
4. Date of admission: ..... Time:.....
5. Lab ID#: .....
6. Specimen type: .....
7. Date of specimen collection: .....
8. Diagnosis: .....

ANNEX 3

AMR-CRF TIER II

ID#: AHC-1709-000

**Hospital/Laboratory Request Form**

AHC  BTB  TKO  SCH

Hospital Service/ ward:.....

Hospital/Laboratory contact: .....

1. Age:..... BoD: .....Days .....Months .....Years
2. Sex:  Male  Female
3. Patient's hospital ID#:.....
4. Date of admission: ..... Time:.....
5. Lab ID#: .....
6. Specimen type: .....
7. Date of specimen collection: .....
8. Diagnosis: .....

## **SPECIMEN COLLECTION – Blood culture (ANNEX 4a)**

### **1. Objective**

To provide instructions on how to perform aseptic blood culture collection ensuring appropriate volume is collected and avoiding contamination of specimen.

### **2. Responsibility**

All phlebotomists (may include nurses, doctors, laboratory staff).

### **3. Principle**

Standard precautions should be observed when drawing blood. The phlebotomist should perform appropriate hand hygiene procedure both pre and post blood collection to avoid possible cross contamination.

### **4. Phlebotomy Procedure (for blood culture)**

- Staff collecting blood will explain the process of blood collection and ask consent from the patient to perform blood collection.
- Staff will prepare the material for specimen collection.
- Label blood culture bottle with patient name, gender, ID number, patient date of birth, site of collection (e.g. left arm), date and time of collection.
- Staff perform hand hygiene
- Staff position the patient
- Put on tourniquet then palpate and select a vein that feels easy to collect blood from. Veins that can be felt are more reliable than those that are seen but not felt well.
- Disinfect the skin using 70% alcohol wipe. Use a circular motion to clean the site first then move outwards to skin beyond the proposed puncture site. If the alcohol wipe is dirty, repeat the procedure until the skin is clean. Discard the alcohol wipe and leave to dry for 1 minute. Clean the area again with betadine 10% (povidone iodine) and leave to dry for 2 minutes.
- Then remove the lid from blood culture bottle and disinfect the rubber stopper using povidone iodine. Leave to dry. It is important to do this prior to organizing needle and syringe to allow time for the povidone to dry (2 minutes).

- Carefully organize the sterile needle and syringe. If syringe comes with incorrect needle size, remove and replace with an appropriate sized needle:
  - Children ( $\leq 14$  years old): using a 23G needle and a 5ml syringe, collect 2-5 ml of blood for a 50ml blood culture bottle
  - Adult ( $> 14$  years old): using a 21G needle and a 10ml syringe collect 10 ml of blood from two different sites for a 100ml blood culture bottle. If the doctor has ordered other tests, use a larger syringe and collect up to 20 mls.
- Without touching the proposed puncture site, carefully insert needle through the skin and into the vein and collect the required amount. If the vein is not entered, you must not use the same needle to re puncture the skin. Replace with a sterile needle.
- Inject total required blood into blood culture bottle and gently invert and mix. If blood is also required for other tests, place blood into tubes only **after** placing in blood culture bottles. Discard the needle and syringe immediately into the sharps container.
- For adults, blood should be collected from 2 sites – 10 mls from each site so repeat the procedure after finding an appropriate vein. Do not take blood from an IV catheter and avoid taking blood proximal to an IV line (may dilute the blood – take blood from another site or distal to the IV catheter). It may be necessary to take blood from the same arm or a leg vein.
- Perform hand hygiene
- Send the blood culture bottles to laboratory along with the microbiology request form that has been completed (signed, date, time). Do not refrigerate.

**Important note for blood culture collection:**

- If possible, collect blood cultures before antibiotic treatment is given to the patient, however do not delay treatment of severely ill patient. If the patient is already receiving antibiotics, take blood culture specimen immediately before the next prescribed dose of antibiotic.
- If the patient has an IV running, draw blood from the opposite arm, other site (leg) or distal to the IV catheter site should a suitable vein be found.
- Blood culture specimen should be requested by the doctors according to the indications agreed to by the hospital.

- The clinician must complete all sections of the request form including clinical syndrome, relevant co-morbidities and whether patient is already on antibiotics. The phlebotomist must sign, date and include time of collection on the request form

**Blood culture bottles**

- Blood culture bottles have either 50ml (paediatric) or 100 ml (adult) broth

## **SPECIMEN COLLECTION – CSF culture (ANNEX 4b)**

### **1. Objective**

To provide instructions on how to perform aseptic CSF culture collection ensuring appropriate volume is collected and avoiding contamination of specimen.

### **2. Responsibility**

The collection of CSF is an invasive procedure and should only be performed by experienced personnel under aseptic conditions.

### **3. Principle**

It is important to adhere to proper biosafety guidelines while handling potentially infectious clinical specimens in order to maintain a safe working environment for patients, health care workers, and laboratorians. Infection may be transmitted from patient to staff and from staff to patient during the procedures described below.

### **4. Preparing for Lumbar puncture**

If possible, three tubes (1 ml each) of CSF should be collected for microbiology, chemistry, and cytology. If only one tube of CSF is available, it should be given to the microbiology laboratory. Because the presence of blood can affect cultures of CSF, if more than one tube of CSF is collected from a patient, the first tube collected (which could contain contaminating blood from the lumbar puncture) should not be the tube sent to the microbiology laboratory.

1. Skin disinfectant: 70% alcohol swab and povidone-iodine
  - Alcohol with concentrations greater than 70% should not be used because the increased concentrations result in decreased bactericidal activity. Do not use alcohol with glycerol added to it.
2. Sterile gloves
  - Be sure to check the expiration date.
3. Sterile gauze
4. Surgical mask
5. Adhesive bandage
6. Lumbar puncture needle
  - 22 gauge/89 mm for adults
  - 23 gauge/64 mm for children
7. Sterile screw-cap tubes
8. Syringe and needle

9. Transport container

**5. Lumbar puncture Procedure (for blood culture)**

Follow all appropriate biosafety precautions.

1. Gather all materials from the CSF collection kit and a puncture-resistant autoclavable container for used needles.
2. Wear surgical mask and sterile latex or nitrile gloves that are impermeable to liquids and change gloves between every patient.
3. Label the collection tubes with appropriate information: patient's name, date and time of specimen collection, and Unique Identification Number. Be sure this number matches the number on both the request and report forms.
4. Ensure that the patient is kept motionless during the lumbar puncture procedure, either sitting up or lying on the side, with his or her back arched forward so that the head almost touches the knees in order to separate the lumbar vertebrae during the procedure (Figure 2).
5. Disinfect the skin along a line drawn between the crests of the two ilia with 70% alcohol and povidone-iodine to clean the surface and remove debris and oils. Allow to dry completely.
6. Position the spinal needle between the 2 vertebral spines at the L4-L5 level and introduce into the skin with the bevel of the needle facing up.  
Accurate placement of the needle is rewarded by a flow of fluid, which normally is clear and colorless.
7. Remove CSF (1 ml minimum, 3-4 ml if possible) and collect into sterile screw-cap tubes. If 3-4 ml CSF is available, use 3 separate tubes and place approximately 1ml into each tube.
8. Withdraw the needle and cover the insertion site with an adhesive bandage. Discard the needle in a puncture-resistant, autoclavable discard container.
9. Remove mask and gloves and discard in an autoclavable container.
10. Wash hands with antibacterial soap and water immediately after removing gloves.

## **AMR Surveillance System Bacterial Isolate Packaging and Transport SOP (ANNEX 5)**

### **1. Objective**

This document describes how to pack, label and prepare documents for transport of bacterial isolates to laboratories within Cambodia for the Cambodia Communicable Disease Control (CCDC) AMR Surveillance System. This ensures isolates are packed and transported in appropriate packaging using the triple packaging systems. It will provide the highest level of safety for transportation.

### **2. Responsibility**

Laboratory personnel

Transportation personnel

### **3. Principle**

The key principles of packaging infectious substances are:

1. To protect the environment
2. To protect the handlers and couriers
3. To protect the sample

#### **Classification of infectious substances**

Dangerous goods are classified and assigned UN numbers.

Infectious substances are classified in Division 6.2.

**UN2814:** Category A infectious substance affecting humans: an infectious substance in a form capable of causing permanent disability or life-threatening or fatal disease in otherwise healthy humans or animals when exposure occurs

Examples: *Burkholderia pseudomallei*

**UN2900:** Category A infectious substance affecting animals

**UN3373:** Category B infectious substance not included in Category A

Examples: *Salmonella* sp., *Staphylococcus aureus*, *Escherichia coli*, *Klebsiella pneumoniae*, *Acinetobacter*, *Haemophilus influenzae*, *Streptococcus pneumoniae*, *Streptococcus agalactiae* (Group B) *Neisseria meningitidis*, *Listeria monocytogenes*, *Streptococcus suis*.

There are no national transport regulations in Cambodia.

Most pathogens transported under the AMR Surveillance System are classified as-

**UN3373:** Category B infectious substance not included in Category A

#### 4. Materials

1. Primary container sealed with parafilm.

Examples: Plastic screw capped TSA tubes (glass tubes should NOT be used) OR agar plates.

2. Leak-proof secondary container

Examples: Screw cap container used for Microbiology EQA in Cambodia. See Figure 1. Right: is a strong, leak proof screw capped container.

**\*\*Attention:** Several primary containers may be placed in one secondary container.

**Figure 1.** Strong leak proof screw capped Secondary container

3. Strong outer package. A large container can hold several Secondary containers if required.

Examples: Large plastic ice box that can be secured with a padlock See Figure 2.

**Figure 2.** Large plastic icebox Strong leak proof screw capped Secondary container

4. Absorbent material to be placed between primary and secondary container.  
Examples: Paper towel or cotton wool.
5. Labels: Hazard and or handling labels
6. Documents: Specimen Shipping List and Referral Request Form

**5. Reagent**

None

**6. Standard and control**

Not required

**7. Sample**

Bacterial Isolates

**8. Procedure**

Isolates are packaged by using **the** triple packaging system.

**8.1.1.** Isolates are inoculated into a media vial or on agar plates

Make sure that the vial or plate is secure by wrapping with parafilm.

This is the primary receptacle.

**8.1.2.** Put the isolate vial or plate in a secure screw cap container (secondary container).

Put absorbent material on the bottom of this container to absorb all fluid in case of breakage.

A Specimen shipping list (F-MCU-045): Is a contents list that records all isolates present in the package. This is placed between the second and third container.

**8.1.3.** Put the secure secondary container in a strong outer package.

**9. Reporting results:**

N/A


**10. Normal Reference Range**

N/A

**11. References**

- IATA infectious substances shipping guideline 10<sup>th</sup> edition, January 2009.
- Guideline on regulations for the transport of infectious substances 2009-2010, World Health Organization.
- DMDP SOP, version 2, document code 018

**12. Flow chart**

### 13. Safety Precaution


Use required PPE when preparing and packaging isolates: lab coat, gloves.

### 14. Supplementary notes

#### a. Packaging:

To ensure isolate, environment and handler protection the triple packaging system is used. Where materials are not available, use similar materials and ensure triple packaging system is followed.

As far as possible, Pack according to Figure 3. with the packaging containers described under **Materials** (see above).

1. Primary receptacle
2. Secondary receptacle (leakproof)
3. Outer container (w/list of itemized contents)

**Figure 3.** Example of a triple packaging system

Source:

<https://web3.unt.edu/riskman/index.php?section=onlinetraining&group=hazmattransportation&module=11>

#### b. Labeling and Marking:

There are 2 types of labelling

##### 1. Hazard labels

Hazard labels indicate the classification of the infectious substance being transported.

See Figure 4.


Figure 4. Hazard labelling for **UN3373**: Category B infectious substance not included in Category A

## 2. Handling labels

See Figure 5. for example, of labels that should be attached to the outer container according to the type of specimen being transported and special handling requirements. These labels can be printed from ANNEXE 2. and attached to the outer container with tape.


Figure 5. Orientation arrow label

### c. Sender and Receiver details

For all shipments include clearly marked on the outer packaging:

1. Sender name, address and telephone number
2. Receiver name, address and telephone number

### d. Documentation: Every shipment must include the following documents

1. **Specimen shipping list** (Refer to the shipping list in the last page): must be placed between the secondary and outer container and list the complete contents of the shipment.
2. **Referral request form**: Clearly record detail of the AMR Surveillance bacterial isolate. One form should be completed for each isolate

### e. Communication: Clear communication between sender, transporter and receiver is essential.

1. The **sender** will always inform the receiver beforehand by **phone** about the moment the shipment is leaving his laboratory and the expected date and time of arrival at the receiver's laboratory. It should be confirmed that there is informed staff available for receiving the shipment.
2. The **receiver** will always confirm the arrival of the shipment to the sender by e-mail or phone.
3. All communication should **be marked clearly** on the Referral request form (F-MCU-044)

<b>NATIONAL AMR SURVEILLANCE</b>	<b>Document code :</b> F-MCU-045	<b>Prepared date:</b> 25/Oct/17	<b>Revision No:00</b>
			<b>Issued date:</b> 28/Oct /17
<b>TITLE:</b> <b>SPECIMEN SHIPPING LIST</b>			<b>Revised date:</b> N/A

**Ship to:** \_\_\_\_\_

**Shipped Date:** \_\_\_\_\_

<b>N°</b>	<b>Patient Name</b>	<b>Patient ID</b>	<b>Specimen Type</b>	<b>Date collection</b>	<b>Test request</b>	<b>Remark</b>
1						
2						
3						
4						
5						
6						
7						
8						
9						
10						
11						
Total						

Prepared by: \_\_\_\_\_

Picked up by: \_\_\_\_\_

Prepared date: \_\_\_\_\_

Picked up date: \_\_\_\_\_

Approved by: \_\_\_\_\_

Received by: \_\_\_\_\_

Approved date: \_\_\_\_\_

Received date: \_\_\_\_\_

**Cambodia-Antimicrobial Resistance Surveillance Form: Isolate Referral Form (ANNEX 6)**

**Patient ID:** \_\_\_\_\_ **Lab ID:** \_\_\_\_\_

**Specimen and Isolate Information**

<input type="checkbox"/> <b>Blood</b> 1. Number of blood culture bottles collected _____ 2. Number of blood culture bottles positive _____ 3. Date specimen collected _____ 4. Date positivity first detected _____ 5. Isolate Identification _____	<input type="checkbox"/> <b>CSF</b> 1. CSF White Blood Cell Count _____ 2. CSF Gram stain _____ 3. Date specimen collected _____ 4. Date culture positive _____ 5. Isolate Identification _____
--	--

6. What method did you use for AST?  Disc Diffusion  Minimal Inhibitory Concentration (MIC)  Other \_\_\_\_\_  
 7. What standard for your AST method was used?  CLSI  EUCAST  CA-SFM Year of standard's last update? \_\_\_\_\_

**Antibiotic Susceptibility Testing (AST) Results** mm-zone diameter  $\mu$ g/mL-MIC  
R-Resistant, I-Intermediate, S-Susceptible: Interpretation

Antibiotics	Disc code	$\mu$ g	mm or $\mu$ g/mL	R/I/S	Antibiotics	Disc code	$\mu$ g	mm or $\mu$ g/mL	R/I/S
Ampicillin	AM	10			Ciprofloxacin	CIP	5		
Amoxicillin/ Clavulanic acid	AMC	20/10			Perfloxacin ( <i>Salmonella</i> ) Surrogate for fluoroquinolones e.g. CIP	PF	5		
Cefazolin	CZ	30			Azithromycin ( <i>Salmonella Typhi</i> )	AZM	15		
Ceftriaxone	CRO	30			Trimethoprim/ Sulfamethoxazole	SXT	1.25/ 23.75		
Cefotaxime #	CTX	30			Penicillin	P	10		
Ceftazidime #	CAZ	30			Cefoxitin ( <i>Staphylococcus</i> ) Surrogate for Oxacillin	FOX	30		
Ceftazidime + Clavulanate #	CAZ/CLA	30/10			Oxacillin ( <i>Staphylococcus</i> )				
Cefotaxime + Clavulanate #	CTX/CLA	30/10			Vancomycin MIC ( <i>Staphylococcus</i> )	V			
Cefepime	FEP	30			Tetracycline	TE	30		
Imipenem	IMP	10			Erythromycin	E	15		
Meropenem	MEM	10			Clindamycin	CC	2		
Gentamicin	GM	10			Penicillin MIC ( <i>Streptococcus</i> )	P			
Amikacin	AN	30			Ceftriaxone MIC ( <i>Streptococcus</i> )	CRO			
Chloramphenicol	C	30			Vancomycin ( <i>Streptococcus</i> )	VA	30		
Other					Other				

**Salmonella serotyping:**  Polyvalent O Positive  Other \_\_\_\_\_

**ESBL** detected?  YES  NO, ESBL Detection Method:  'Keyhole' effect  Clavulanate combination disc testing  
 If Clavulanate Combination disc testing was performed, please ensure zone sizes are recorded for # antibiotics  
**ICR** (Inducible Clindamycin Resistance) detected?  YES (D-test Positive)  NO (D-test Negative)  
 Is isolate Resistant or Intermediate to Imipenem or Meropenem?  YES  NO  
 Did you test for **Carbapenamase**?  YES  NO  
 Detection method:  mCIM (modified Carbapenam Inactivation Method)  Other .....What was the result?.....

Completed by: \_\_\_\_\_ Date of Referral: \_\_\_\_\_ Contact Number: \_\_\_\_\_

**ANNEX 7**

**KERRY EXPRESSS, CAMBODIA**

Head Office: # 38, Street 86T (Near Wat Samrong Andeth)

Sangkat Phnom Penh Thmei, Khan Sen Sok,

Phnom Penh, Cambodia

Tel: 023 231 232 (Receptionist)

Tel: 081 555 199 (Sale service)

**EMS – EXPRESS MAIL SERVICE, CAMBODIA**

No. 486BC, Street. Mao Tse Toung Blvd (245)

Phnom Penh, Cambodia

Tel: +855 23 88 87 74

Mobile: +855 88 87 11 839

Annex 8